



## IN THIS ISSUE


- Oticon Amigo™ Personal Communication System
- Sennheiser™ Set 810
- Printacall Home Room Loop System
- Oricom TP1100 Volume/Tone Control Phone
- Bellman™ Audio Maxi communicator
- Alarm Clocks
- Bellman™ flashing & vibrating Smoke Alarm
- TTY communication options from Telstra

Call  
for our  
latest  
brochure

“Improving the life  
of Deaf and hearing  
impaired people  
since 1984.”

## Introducing the Oticon Amigo™ System -A fully portable FM, go anywhere companion.

With the Oticon Amigo™ your direct access to sound without annoying background noise will bring you all the clarity you desire. The NEW Oticon Amigo™ supersedes the previously well known personal Easy Listener FM system which reputedly delivered quality sound with a multitude of applications. This new system has all the original qualities of the Easy Listener with the added bonus of being more compact in size and features a multi channel facility.



### When would I benefit by an Amigo™ System?

If you struggle in any of the following situations then you should seriously consider an Oticon Amigo™ System.

#### Life is too short to miss out on anything.

- Do you strain to hear in a classroom / lecture / church situation when you are some distance from the speaker and sound has to travel?
- When travelling in a car do you experience difficulty communicating?
- Is it difficult to hear when attending functions / club meetings with an around the table venue?
- Do you struggle to hear people at family /group gatherings?
- Is it difficult to hear your partner when you dine out?

The Amigo™ is compatible with all the optional accessories that were previously available for the Easy Listener:

- A neckloop for the hearing aid T switch user
- Headset, stetoclip or earbuds for the non hearing aid user
- A conference microphone for around the table meetings

### Interested and want to know more?

For further information on which Oticon Amigo™ System would best meet your individual situation speak to one of our friendly, knowledgeable staff in our showroom, by phone or email. They will be pleased to discuss your needs and advise you.

\* Trademark of Oticon A/S.  
Phonic Ear / Logia Denmark

# Did You Know About?

## Captioned TV - The TV listening solution?

In order to view TV with captions you can use either a teletext TV, or if you live in an area which receives digital broadcasts you can either use a digital TV or digital set-top box. Both of these digital devices should enable you to view captioned programs. Be aware that not all set top boxes can access captioning so we recommend you ask for a demonstration before buying.

Watching TV doesn't need to be hard work. Although Printacall doesn't sell TVs or captioning equipment they

do have other TV listening solutions. Ask Printacall about their range of cordless Sennheiser products (that do not disturb others) or how TV can be enjoyed without having to wear equipment (with the installation of an audio induction loop (room loop) in your home).



## Improve your TV Listening

### A Printacall Home Room Loop System brings the sound to you


A home room loop enables a person wearing a hearing aid with a T switch to access the sound from a TV, radio or other equipment via the loop amplifier. The loop minimises distracting background noise from within the room and from outside.

The volume of sound from the TV or radio through the room loop to the user is at a more natural level. The whole family can enjoy the sounds of entertainment without uncomfortable noise levels.

The user of the room loop is free to move anywhere within the looped room and enjoy the sounds at their individual volume level, without interfering with other people.


#### PRINTACALL LOOP SYSTEMS

Does your church or meeting hall have an audio loop system? Printacall is a specialist in the design and installation of powered audio induction loop systems with over 20 years of experience. Our skilled technicians have installed hundreds of loop systems into public access buildings.



### Set 810 Mono/ Stereo Infra-Red System

Turn up the volume when watching TV without disturbing others. Connect the transmitter to the audio output socket on your TV, then adjust the volume on the headset.



Also available in a neck loop style to deliver high-level quality sound from the TV or radio to your hearing aid with a "T" switch.

\* Trademark of Sennheiser electronic GmbH & Co KG.

## Oricom TP1100™ Loud and Clear Amplified Telephone

The TP 1100 is the ideal telephone solution for those with mild-to-moderate hearing loss and for use in noisy environments. It's adjustable

tone and volume control, telecoil compatibility and 12 memory keys makes conversing easy.

\* Trademark of Oricom International Pty Ltd.



## ORDERING OPTIONS

BY PHONE  
02 9809 2392

BY TTY  
02 9809 1283

BY FAX  
02 9809 2345

BY EMAIL  
sales@printacall.com.au

BY POST  
2 Doig Ave., Denistone  
East, NSW, 2112

For more information visit  
www.printacall.com.au

# Products

## Bellman™ Audio Maxi Personal Communicator

The perfect solution for hospitalised hearing impaired individuals who, for whatever reason, are unable to utilise hearing aids.

The user simply slips on the headset or stetoclip and aims the Maxi at the speaker to receive quality sound. Tone and volume controls on the Maxi individualise its performance and battery operation makes it very transportable.

A recent happy recipient of this

equipment can't speak highly enough of it:

*"It has been the most wonderful thing for my husband he would have been battling without it. I wish the hospital would buy some. During his stay he has been bombarded with questions needing answers from not only the doctors and nurses, but physiotherapists, radiologists, menu and cleaning staff to name a few. This equipment enabled him to reply and literally put him on the air. This has been phenomenal. Without it he would have been up the*



A neckloop option is also available for hearing aid "T" switch users.

*proverbial creek! The nurses are saying they wish the hospital had access to this equipment as it would make their lives much easier when nursing elderly hearing-impaired patients."*

I. Watson

TM Trademark of Bellman + Symfon AB, Sweden.

## Shake Awake™ Alarm Clock

Shake Awake™ is a compact, self contained, battery powered vibrating alarm clock. No wiring or separate vibrator is required. It nestles under the pillow or in the pocket.

This advanced design affords the hearing impaired person mobility and convenience. It can be used as a wake up device or as a cooking or other event timer.



TM - Trademark of FOB Inc.

## Vibralarm™ Clock with Bed Shaker

This digital clock lets you set the loudness and the tone of the alarm that best suits your needs. Attach a bedshaker (sold separately) to place under your pillow or mattress to shake you awake. Heavy sleeper? Use both the bedshaker and the audible alarm together. This is the alarm clock you need when you don't want to oversleep.



### FEATURES

- Extra loud alarm clock.
- Audible alarm with variable volume and tone.
- Auxiliary jack to attach Bed Shaker (sold separately).
- Night light and snooze control.
- Selectable controls for steady or pulse vibration, with or without audible alarm.
- Battery back up (9V battery not included).

TM - Trademark of Global Assistive Devices.

## Wake 'N' Shake Alarm Clock

A jumbo alarm clock and shaker ideal for people with a hearing impairment or for those who don't want to wake their partner in the morning. It can also be connected to a phone to provide a loud ringer, flashing strobe light and vibrating shaker, making it easier to hear and see when you have an incoming call or an alarm alert.

Wake 'N' Shake has an oversized display for easy reading.



### FEATURES

- Large easy to read display.
- Large vibrating alert pad.
- Extra bright flashing strobe.
- Select vibrate/alarm/strobe.
- Extra loud adjustable volume (+95dB max).
- Tone control for alarm.
- Connects to any phone and detects ring.
- Simple to connect.

TM - Trademark of Oricom International Pty Ltd.

# The Bellman™ Visit Smoke Alarm

The Bellman™ Visit system helps people in difficult listening situations. The Smoke Alarm transmitter is activated by the presence of smoke in the event of a fire, transmitting a signal to the various Bellman™ Visit receivers, such as the Visit Flash, which alerts the user by light, and the Visit Bedshaker, which alerts by vibrating.



Visit Bedshaker

Visit Smoke Alarm

## FEATURES:

- **Easy to install** - simply connect the flash to a power point and fit one or more smoke alarms to the ceiling.
- **Portable** - battery backup allows for wire less operation of the Flash for up to 3 days.
- **Dual Signals** - provides an alert with a strobe light and vibrating pad for night safety. The Flash and Bed Shaker will respond as long as the Smoke Alarm is activated.
- **Good Range** - range between the Flash and the Smoke Alarm may be up to 30 metres subject to structural conditions. Thick walls will decrease the range.
- **Easy Testing** - a test button on the alarm will give you peace of mind and a low battery indicator will alert you when to replace the battery.

TM Trademark of Bellman + Symfon AB, Sweden.

# TTY communication options from Telstra

Teletypewriters (TTYs) allow customers to send and receive text conversation over the telephone network. TTYs provide the flexibility to conduct text conversations with other TTY users, or to use the National Relay Service to conduct text to voice calls to other callers.

Telstra offers eligible customers the option of renting one of three TTY models under its extensive Disability Equipment Program.

For further information please contact:

- Telstra's Disability Enquiry Hotline
- FREECALL\* 1800 808 981 (TTY),
- FREECALL\* 1800 068 424 (Voice),
- or Email: Disability@online.telstra.com.au

\*A FREECALL from most fixed phones. Printacall is contracted to Telstra to supply, distribute and maintain TTYs.

## Uniphone 1150™

The telephone for everyone, combination TTY, telephone and amplified phone all in one!

\* Trademark of Ultratec Inc., Madison, WI USA



## FEATURES:

- Turbo Code™ with interrupt facility
- Voice carry-over (VCO) and Hearing carry-over (HCO)
- Built-in ring flasher for incoming calls
- Hearing-Aid coupler
- Volume control on handset
- Voice announcer. ASCII
- 48 character LCD display

PRODUCT SPECIFICATIONS AND AVAILABILITY MAY CHANGE WITHOUT PRIOR NOTICE

## Privacy Statement

The personal information which we collect includes names and addresses and, in the case of customers, details of products and services supplied. Such information is recorded in one or more databases of:

- Our customers who have bought equipment or services, some of whom have lodged warranty cards with us.

- Persons who have asked to go onto our mailing list for product information and newsletters.
- Businesses, organisations and government bodies, and individuals connected with them, whose details have been compiled from telephone directories or subsequent telephone contact.

We regard our database contents as confidential to Printacall

and our staff is instructed to preserve such contents as confidential.

We do not disclose personal information to any outside third party unless:

- it is contracted to us to produce printing and mailing services, and in such cases we require the third party to be bound by similar privacy rules to our own;

- required by law to do so; or
- consent is given by the person about whom the information is held.

You may at any time, on request, have access to the personal information we hold about you and you may ask us to correct or update it. Please let us know if you do not wish to receive marketing information or if you would like a copy of our Privacy Policy.



PH: 02 9809 2392

- TTY: 02 9809 1283
- EMAIL: sales@printacall.com.au
- FAX: 02 9809 2345
- WEBSITE: www.printacall.com.au
- 2 Doig Ave.,
- Denistone East,
- NSW, 2112