

T-Link™ for Mobile Phones - Hands Free Kit for Hearing Aid Users

A must for Hearing Aid Users

FEATURES

- Allows for interference-free access to most Mobile Phones.
- No batteries or user serviceable items required.
- Both Binaural and Monaural T-Links™ available.
- Simply place the hook behind the ear fitted with a hearing aid and switch to the "T" position.



TM - Trademark of Sensorcom, UK.

Call **PRINTACALL** for Service

PRODUCT SPECIFICATIONS AND AVAILABILITY MAY CHANGE WITHOUT PRIOR NOTICE

A division of All Mod Comms Pty Ltd

ACN 002 863 103

PRINTACALL



2 DOIG AVE., DENISTONE EAST, NSW 2112

email: sales@printacall.com.au

web: www.printacall.com.au

TEL : (02) 9809 2392

FAX : (02) 9809 2345

TTY : (02) 9809 1283